

CT Medical Images Dataset Analysis

A Detailed Report on Pixel Intensity Analysis

By

**Gopalakrishnan Kumar, IIT-B Alumnus,
Freelance Data Science Consultant.**

LinkedIn URL –

**[https://www.linkedin.com/in/gopalakrishnan
-kumar-a73301110/](https://www.linkedin.com/in/gopalakrishnan-kumar-a73301110/)**

Git Repository

**[https://github.com/Gopalakrishnan-
Kumar/Python-for-Data-Science](https://github.com/Gopalakrishnan-Kumar/Python-for-Data-Science)**

Website/blog URL

<https://www.kaggle.com/gopalkk1>

Introduction to CT Medical Imaging

- Computed Tomography (CT) scans are a cornerstone in modern medical diagnostics, providing detailed cross-sectional images of the body. These images are created by combining multiple X-ray measurements taken from different angles, allowing physicians to examine internal organs, bones, soft tissues, and blood vessels with high precision. The ability to visualize internal structures in detail makes CT scans particularly valuable for diagnosing a wide range of conditions, including cancers, cardiovascular diseases, and musculoskeletal disorders.

Introduction - Applications in Machine Learning

- CT medical images provide a rich source of data for developing machine learning models that can assist in medical diagnostics. Some key applications include:

Introduction - Classification

- Machine learning models can be trained to classify CT images based on the presence or absence of specific medical conditions. For example, a model could learn to distinguish between healthy and cancerous tissues.

Introduction - Segmentation

- Segmentation involves partitioning the image into regions corresponding to different anatomical structures or abnormalities. This is particularly useful for isolating tumors, lesions, or other areas of interest.

Introduction - Diagnosis

- AI-driven diagnostic tools can analyze CT images to provide diagnostic suggestions, helping radiologists make more informed decisions. These tools can highlight areas of concern, suggest potential diagnoses, and prioritize cases that require urgent attention.

Introduction - Tumor Detection

- Detecting and classifying tumors in the brain using CT images is a critical application of this technology. Early and accurate detection can significantly improve patient outcomes.

Objectives of the Study

- To study the Pixel Intensity analysis of the given CT Scan Image in terms of
 - a. Normalization
 - b. Noise Reduction
 - c. Contrast Enhancement
- To enable medical professionals to identify the type of tissues in the CT Scan image using pixel intensity and detect abnormalities

Dataset Overview


- The CT Medical Images dataset from Kaggle is a comprehensive collection of CT scan images in TIFF format, accompanied by metadata that includes patient age, contrast usage, and other relevant details. This dataset is designed for use in training machine learning models aimed at detecting and diagnosing medical conditions based on CT scan data. By leveraging such datasets, researchers and developers can create algorithms that assist radiologists in making more accurate and timely diagnoses.

Sample Image Details

- One of the files in this dataset is 'ID_0000_AGE_0060_CONTRAST_1_CT.tiff'. This file contains a CT scan image of a 60-year-old patient, taken with contrast enhancement. The metadata for this image is as follows:
 - **ID:** 0000
 - **Age:** 60
 - **Contrast:** Yes (1)
 - **Modality:** CT
 - Contrast enhancement in CT scans involves the use of contrast agents to improve the visibility of blood vessels, tissues, and other structures, making it easier to identify abnormalities.
 - **Location:** Chest, **Species:** Human, **Condition:** Lung Cancer
 - **Number of the Patient:** 1

Lung Cancer Image

TIFF Image



Technical Analysis

- CT images in this dataset are high-resolution TIFF files, providing detailed visual information necessary for medical analysis. The specific image 'ID_0000_AGE_0060_CONTRAST_1_CT.tiff' has the following technical specifications:
- **Resolution:** High-resolution TIFF image
- **Image Dimensions:** Typically, 512x512 pixels (the exact dimensions should be confirmed from the file properties)

Methodology


Data Preparation

- The first step involves preprocessing the images to ensure consistency and enhance the relevant features. This includes normalization, noise reduction, and contrast enhancement. Data augmentation techniques, such as rotating and flipping the images, can also be used to increase the diversity of the training dataset.


Preprocessing Steps

- Before using the images for machine learning, several preprocessing steps are often necessary:
- **Normalization:** Adjusting the pixel intensity values to a common scale.
- **Noise Reduction:** Applying filters to remove any noise that might obscure important details.
- **Contrast Enhancement:** Further improving the visibility of specific structures.


Normalization



Noise Reduction



Contrast Enhancement




Visual Analysis

- The image 'ID_0000_AGE_0060_CONTRAST_1_CT.tiff' provides a clear visualization of the patient's internal structures. Upon examining the image, one can observe various anatomical features, including the enhanced blood vessels due to the contrast agent. These areas of higher contrast are crucial for identifying potential issues such as blockages, tumors, or other abnormalities.
- In this image, regions with increased contrast appear brighter, highlighting the blood vessels and other structures that have absorbed the contrast agent. This contrast differentiation is essential for distinguishing between different types of tissues and identifying anomalies.


Data Analysis- Pixel Intensity

- The pixel intensity values in the CT scan represent different tissue densities. Higher intensity values correspond to denser tissues such as bones, while lower intensity values represent softer tissues. By analyzing the distribution and range of pixel intensities, medical professionals can differentiate between various tissue types and detect abnormalities.


Boxplot of Pixel Intensities




Scatterplot of Pixel Intensities



Bar Chart of Pixel Intensities





Model Training

The "sequential_4" model represents a convolutional neural network (CNN) designed for image classification tasks. With a total of 14,835,234 trainable parameters, this model follows a standard architecture for processing and classifying images. It begins with convolutional layers to extract features from input images, followed by max-pooling layers to reduce spatial dimensions. The model then flattens the feature maps into a one-dimensional array before passing them through dense layers for classification. The final output layer employs a softmax activation function, suitable for multi-class classification tasks. Overall, this CNN architecture is structured to learn and classify complex patterns in medical images, making it suitable for tasks such as identifying and diagnosing medical conditions based on CT scan images.

Model: "Sequential_4"

Layer (type)	Output Shape	Param #
=====		
conv2d_9 (Conv2D)	(None, 254, 254, 32)	896
max_pooling2d_9	(MaxPooling2 (None, 127, 127, 32)	0
conv2d_10 (Conv2D)	(None, 125, 125, 64)	18496
max_pooling2d_10	(MaxPooling (None, 62, 62, 64)	0
conv2d_11 (Conv2D)	(None, 60, 60, 128)	73856
max_pooling2d_11	(MaxPooling (None, 30, 30, 128)	0
flatten_3 (Flatten)	(None, 115200)	0
dense_6 (Dense)	(None, 128)	14745728
dense_7 (Dense)	(None, 2)	258
=====		

Total params: 14,835,234

Trainable params: 14,835,234

Non-trainable params: 0

1. Input Layer (Conv2D):

1. Input Shape: (image_height, image_width, 3)
2. Output Shape: (None, 254, 254, 32)
3. Parameters: 896
4. Activation Function: ReLU

2. MaxPooling2D Layer:

1. Pooling Size: (2, 2)
2. Output Shape: (None, 127, 127, 32)

3. Conv2D Layer:

1. Output Shape: (None, 125, 125, 64)
2. Parameters: 18,496
3. Activation Function: ReLU

4. MaxPooling2D Layer:

1. Pooling Size: (2, 2)
2. Output Shape: (None, 62, 62, 64)

5. Conv2D Layer:

1. Output Shape: (None, 60, 60, 128)
2. Parameters: 73,856
3. Activation Function: ReLU

6. MaxPooling2D Layer:

1. Pooling Size: (2, 2)
2. Output Shape: (None, 30, 30, 128)

7. Flatten Layer:

1. Output Shape: (None, 115200)

8. Dense Layer:

1. Output Shape: (None, 128)
2. Parameters: 14,745,728
3. Activation Function: ReLU

9. Output Layer (Dense):

1. Output Shape: (None, 2)
2. Parameters: 258
3. Activation Function: Softmax
(assuming 2 output classes for binary classification)

The total number of trainable parameters in the model is 14,835,234. This model follows a typical **convolutional neural network** architecture with convolutional layers followed by max-pooling layers, flattening, and dense layers for classification.

Conclusion and Future Work

- The Location was in the Chest. The condition was Lung cancer. The analysis of CT medical images using machine learning holds great promise for improving diagnostic accuracy and efficiency. By leveraging high-resolution CT scans and advanced AI algorithms, it is possible to assist radiologists in making more accurate and timely diagnoses. Future work in this area could involve exploring more advanced models, integrating CT images with other imaging modalities (such as MRI or PET scans), and collaborating with medical professionals to ensure the models are clinically relevant and effective.

Reference

Kaggle